

# USER MANUAL



Pls read instructions carefully  
before operation



2 in 1  
Vacuum Cupping Massage

## Contents

1. Safety Notice .....	2
2. Vacuum Breast Enlargement Cupping .....	3
3. Machine and Accessories.....	4
4. How to Operate? .....	8
5. Device Daily Maintenance .....	11
6. Technical Parameter & Packing List.....	11
7. After-sale Service & Maintenance.....	12
8. Warranty Card .....	13

## 1. Safety Notice

1. Please don't disassemble the instrument which hasn't been indicated in our instruction. All the repairs should be done by authorized professionals.
2. Please don't assemble or operate the instrument near wet place. Water and any liquid are not allowed to pour on screen.
3. Please pull off the plug when you find conditions below.

The equipment touched liquid:

- A. The equipment has abnormal smell, smoking or loud and strong noise.
  - B. Cable is broken
  - C. The equipment is fallen down and broken
4. No stuff on the cable, keep people from stepping.
  5. Please pull off the power cord after finish operation.
  6. Don't put anything into the equipment's jack, otherwise it may cause fire or touch the electricity. If something enters the equipment, please don't pick it out by yourself, contract distributor or our company immediately
  7. Please don't put the equipment on the anything that can't stand steadily, avoid it turns over and damage.
  8. Our company has the right to modify the information in this introduction before notify! And keep the explanations right for the above.

## 2. Vacuum Breast Enlargement Cupping

This instrument is the most effective and direct breast building instrument, which develops under the basis principle of physiology and physical sucking theory. And its slight but concentrated effect could exercise the muscular tissue mostly (breast muscle develops----next to neural tissue----finally reach hypophysis) and stimulate hypophysis to excrete hormone, thus to exercise muscular tissue stretchy as well as fibre function, also enlarging breast naturally and effectively. It suits Eastern feminine physiology and proved by years of clinical confirmation that it does no harm to human body but does effect for breast beauty.

### Function

**Breast care:** Breast plumping & lifting、Correct nipple & chest、Recover breast elasticity.

**Facial care:** Face thinning、Remove nose blackhead、acne、eye circle、 eye bag、 blood marks、 double chin、 & Lifting canthus.

**Back care:** Negative pressure metabolism、 Remove measles、 Press cavity、 Toxin excrete、 Tighten skin、 Improve skin elasticity、 Lymphs toxin removal.

**Buttock care:** Butt lifting, butt firming


**Fatness care:** Losing weight、 Slimming、 Tightening abdomens、 Slimming waist.

### 3. Machine and Accessories





Total 32pcs cups including 4 pairs breast cups.




Air Release Lid

Vibration On/Off

DC Line

# KLD Skin Care



12 tubes and connectors

Two short tubes

DC line for red light and vibration on breast cups

Tee regulator air control

Vacuum heating handle

Main long tube





## 4. How to Operate?

1. Clean the skin.
2. Smear slimming oil or tightening oil on the skin.
3. Use Vacuum heating handle to massage body for 5-15 minutes.
4. Use proper cups to treat breast, back, leg, arm, abdominal.

Adjust proper suction time for customers.


Usually treatment time 20-30 minutes.

After treatment, redness is normal.

When massage face, do not use big suction time, reduce the suction time, avoid redness.

Please stop the machine after one hour working, let it rest for 30 minutes.

About installation and operation details, check our youtube video.



**Massage breast before treatment**






Facial massage directions



Vacuum heating handle (body massage, detoxing)





Scratch  
detox



Cupping  
Comfortable meridian



Cupping  
Relieves fatigue



Cupping  
detox

## 5. Device Daily Maintenance

1. Clear the device with warm and clean water in time after use
2. Do not clean the device with anything contains alcohol.
3. Please power off the device after use.
4. Check the electricity carefully before every operation.
5. Please use the specified accessories provided or recommended by the original manufacturer to ensure the effect and device life.
6. Keep the device away from humid or very hot place.
7. Keep the device away from direct sunlight.
8. Do not drop, strike, beat, throw the device and handle,  
Or else it may cause damage.

## 6. Technical Parameter & Packing List

### Technique Parameters:

<b>Voltage</b>	100V~240V
<b>Power</b>	250W
<b>G.W.</b>	5KGS
<b>Packing Size (cm)</b>	50*35*32CM

### Packing List

Item Name	Unit	No.	Remark
<b>Main Machine</b>	unit	1	
<b>Vacuum heating Handle</b>	pc	1	
<b>Breast cups</b>	pairs	4	
<b>Other cups</b>	pc	24	
<b>DC line</b>	pc	1	
<b>Power cord</b>	pc	1	
<b>Tubes</b>	pc	15	
<b>Air tee regulator</b>	pc	1	

## 7. After-sale Service & Maintenance

The term of maintaining guarantee is main device 1 year and spare part 3 months from the date of purchasing according to the receipt. (Expendable is not included in it.) Others not included in maintaining guarantee or over the time, we also supply to maintain and renew but need charges. Please contact with the dealer if need change of purchasing the parts and easily expendable products.

**A.** We will guarantee the quality of our equipment, normally there is no quality problem

1. If there is a deficiency, please email or call us on the reasons, it would be best to dealt with by email or telephone. Because some of the problems are caused by using the machine in a wrong way.
2. We suggest the client to take photos of the area where the problem occurs, we will pass them to our technician , our technician will give a satisfactory answer within 1-7 working days
3. If the problem comes from the spare part within the warranty, we could consider in changing the spare part to our client, which could be sent with his/her next order. If the next order is uncertain, we could pay freight cost to send the spare parts.
4. If the problem comes from the main device within the warranty, we could advice our client to post the device back to our company for repair by ordinary post, or we could re-send a new one to our client. If the cost of posting the problematic device back to BEST is more than the value of a new device, BEST would take it into consideration, according to the actual status; apply to re-send a completely new device back to our client.

**B.** We won't supply maintaining guarantee if troubles as below:

- a) Disassembling the machine by self;
- b) Fall it without caution;
- c) Lack of reasonable keeping and protect;
- d) Not operate according to the user manual.
- e) Expendable reasons

If changes of the design and specification of the machine, we will not notice you in addition. If any queries, please feel free to contact with the local marketing company.

All in one, we will do our best to serve you and also need your fully cooperation.

## 8. Warranty Card

Product Name:		
Purchase Date:		
Deadline:		
Seller Name:		
Seller Add:		
Seller Tel:		
Repair Date:	Fault Manifestation	Broken Stuff



# Quality & Service Our Focus